



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Streak Retinoscope

Product Code :
NMCBS137029



Description :

Streak Retinoscope is exceptionally bright. It is a Streak and Spot Retinoscope in one and is converted by simply changing the bulb. Retinoscope with D size battery handles in case or zipper pouch packing. Rotation and Vergennes movements independently operated for those who prefer separate controls. Continuous 360° rotation for easy and accurate positioning of streak for axis determination. Streak Retinoscope Streak Retinoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,