



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Xenon Light Source

Product Code :
NMCBS137031



Description :

Xenon Light Source We present luminous xenon light source in the marketplace. We make available different ranges and specification adding choice to your exact requirements. Like our pulsed xenon light source provides crystal clear light, without any blurring of the image. These xenon light source is one of the most important equipments for any medical application, either in the Operation Theater or clinic. Further, the endoscopic xenon light sources we supply in the market are variedly used for endoscopic surgeries as they provide slender but bright light. It really becomes easy to have a clear view and image of the affected part with the help of our xenon light sources. Moreover, we are also one of the leading xenon light sources manufacturers of the country. Xenon Light Source, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)