

**Product Name :**  
Applanation Tonometer**Product Code :**  
NMCBS137036**Description :**

Applanation Tonometer Applanation Tonometer is produced for those who wish the tonometer to remain permanently on the slit lamp. It is mounted on a pivot on the microscope and for examination is swung forward in front of the microscope. The appellation surface is to be examined in a monocular manner through the left ocular. Good service Measured interocular pressure Supply a full series high quality optical/ophthalmic/medical instrument Specifications: Measuring Range: 0~10.6kpa Movement of light circle:  $1.53 \times 2 = 3.06\text{mm}$  Diameter of Probe: 7mm Moving angle of probe: 3mm With Slit Lamps Applanation Tonometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers