

**Product Name :**  
DRYING OVEN (TRAY DRIERS) ELECTRICAL HEATED &  
STEAM HEATED

**Product Code :**  
NMCBS138012



**Description :**

DRYING OVEN (TRAY DRIERS) ELECTRICAL HEATED & STEAM HEATED (Maximum temperature 100 degree centigrade). Generally as per enclosed catalogue. (A) ELECTRICALLY HEATED ; The unit is electrically heated with 'ESCORTS' U tube heaters and provided with thermostatic control dial thermometer. [B] STEAM HEATED : By means of Steam Radiators (M.S. Hot Dipped Galvanized). The unit is without temperature controller, steam fittings valves etc. DRYING OVEN (TRAY DRIERS) ELECTRICAL HEATED & STEAM HEATED, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,