



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
VACUUM TRAY DRYER

Product Code :
NMCBS138013



Description :

VACUUM TRAY DRYER 1/2 Overall body of thick S.S. Plate with optional SS 304Q or SS 316Q. 1/2 Option of hot water circulation by limped coil for outer body to avoid dew/moisture formation on the internal chamber of dryer. 1/2 Option of providing receiver condenser or vacuum pump to recover the harmful chemicals or expensive solvents 1/2 Maximum contact surface for efficient heat transfer to product by giving flattened bottom of the trays & top plate of heating shelves. 1/2 Longer & durable heating pad plates made from thick SS sheet fully argon arc gas Welded & provided with S.S. 304Q square rod stiffeners 1/2 Hinged 'C' Clamp with knobs all across the door to avoid leakages VACUUM TRAY DRYER, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,