



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
VACUUM TRAY DRYER

Product Code :
NMCBS138013



Description :

VACUUM TRAY DRYER 1. Overall body of thick S.S. Plate with optional SS 304Q or SS 316Q. 2. Option of hot water circulation by limped coil for outer body to avoid dew/moisture formation on the internal chamber of dryer. 3. Option of providing receiver condenser or vacuum pump to recover the harmful chemicals or expensive solvents 4. Maximum contact surface for efficient heat transfer to product by giving flattened bottom of the trays & top plate of heating shelves. 5. Longer & durable heating pad plates made from thick SS sheet fully argon arc gas Welded & provided with S.S. 304Q square rod stiffeners 6. Hinged 'C' Clamp with knobs all across the door to avoid leakages VACUUM TRAY DRYER, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE