

**Product Name :**  
STAINLESS STBEL CYLINDRICAL TANKS**Product Code :**  
NMCBS138021**Description :**

STAINLESS STBEL CYLINDRICAL TANKS Made out of stainless steel (304/316 q.) 16 gauge with all the welding done by argon arc process. The tanks are polished and buffed from inside and outside to give mirror finish and are supplied complete with lid, one outlet valve with hose nipple, S.S. cladded M.S. Stand having castor wheels etc. Model T-50 and T-100 are supplied without valves and stand. The dimensions are approximate and are subject to change depending on availability of S.S. Sheets STAINLESS STBEL CYLINDRICAL TANKS, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

**Technical Specification :**

Made out of stainless steel sheets (304 Quality) with welding done by Argon Arc process.

The tanks are polished and buffed from inside and outside to give mirror finish and are supplied complete with lid I outlet valve with nose nipple M.S stand with castor wheels etc.

These tanks are used in industries for storage system.

