

Product Name :
RECATION VESSEL**Product Code :**
NMCBS138022**Description :**

RECATION VESSEL Stainless Steel 304/316L quality. $\frac{1}{2}$ MS/SS Jacket or limpet coil for steam heating / cooling. $\frac{1}{2}$ Reaction Vessel can be provided with different types of Agitator like Anchor type / Propeller type / Paddle type Pitched blade turbine type of with high speed homogenizer. $\frac{1}{2}$ Shaft entry sealed by Oil seal / Plain or Water cooled Stuffing Box / Single or Double Mechanical seal. $\frac{1}{2}$ High density glass wool insulation & S.S. Cladding can be done on demand.. RECATION VESSEL, Pharmacy Equipment exporters, Pharmacy Equipment manufacturers

Technical Specification :

Type: Reaction Vessels
Material: Stainless Steel
Shape: Cylinder Shape, Round
Color: Brown, Golden-Grey, Grey, Light Green, Metallic, Shiny Silver, Silver
Surface Finishing: Chemical Coated, Powder Coated
Feature: Anti Corrosive, Durable, Eco-Friendly, High Quality, Shiny Look, Stable Performance

