

**Product Name :**  
Binocular Loupes**Product Code :**  
NMCBS139009**Description :**

Binocular Loupes We offer a wide array of superlative binocular loupes in various sizes and dimensions. We have carved a niche as one of the salient ophthalmology binocular loupes manufacturers, based in India. Our binocular loupes are manufactured from the genuine quality raw materials so as to ensure their excellent performance. Lightweight precision optics is used in our dental binocular loupes. They find a wide area of application in various fields like dermatology, ophthalmology, general medicine, plastic and neurosurgery, dental medicine, veterinary, ENT etc. Binocular Loupe integral with precision Optics with exceptionally lightweight. They are perfect for many applications. General Medicine, Dermatology, General, Plastic and Neurosurgery, ENT, Ophthalmology, Veterinary and Dental Medicine etc AEI-12.010 Binocular Loupe 3.5x with Small Frame AEI-12.020 Binocular Loupe 3.5x with Large Frame AEI-13.010 Binocular Loupe 4.5x with Small Frame AEI-13.020 Binocular Loupe 4.5x with Large Frame AEI-13.200 Binocular Loupe 2.0x with Frame Small/Large AEI-13.250 Binocular Loupe 2.5x with Frame Small/Large Binocular Loupes, ENT Equipment exporters, ENT Equipment manufacturers