



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Anti Bedsore Mattress

Product Code :
HBMAT2010003



Description :

Anti Bedsore Mattress:
Anti-Bedsore mattress (Air Jet Type).
Tiny Holes on the surface for decreasing humidity.
Replaces the foul air and reduce growth of microbes.
It can effectively prevent from bedsore by dispersing weight of patient's body.
Suitable for long term bedridden patients.

Technical Specification :

Very easy to operate - Adjustable comfort range with pressure-adjustable control knob
Provides perfect support to bedridden patients
Waterproof, antibacterial and fire resistant mattress offers maximum benefits
It relieves the pain of sores and pressure ulcers by maintaining alternating pressure
Perfect choice while you are looking for a portable, light weight and a perfect space saver option



Medical Hospital Equipments Manufacturers

Website: SITE