

Product Name :
Artificial Eye**Product Code :**
NMCBS139019**Description :**

Artificial Eye We are one of the eminent organization engaged in offering quality assured range of Cosmetic Shell (Artificial Eye) which is used in case of disfigured eye. Moreover, it can only be used if the patient do not have any danger in the retention of sightless eye. Made of acrylic, these are lighter and thinner than the artificial eye. Available in universal size of 22mm, 24mm, 26mm and 28mm these are extensively used after evisceration. There are two type of cosmetic shells : Cosmetic shell with transparent sclera: This is used in the absence of iris or lack of iris pigmentation when the sclera portion of the eye is not damaged and only the iris is disfigured. Cosmetic shell with opaque sclera: This is used after evisceration when the sclera portion is damaged along with the iris portion. Specifications: Material: 100% natural silk. Style: Solid silk satin fabric with embroidered logo elastic lace band Size: 18.5 x 8.5 cm Packing: One piece in a poly bag Carton size: 62 x 40 x 50 cm MOQ: 200 pcs / color Sample lead time: 7 days Color: Decided by customers Payment terms: TT , L/C, Paypal or Western Union. Delivery time: 30 days after receipt of the deposit or confirm the order Artificial Eye, ENT Equipment exporters, ENT Equipment manufacturers