

**Product Name :**  
Prism Bar Set**Product Code :**  
NMCBS139020**Description :**

Prism Bar Set Catering to the demands of clients we are engaged in supplying quality range of Prism Bar Set. The whole set includes cross cylinder, retinoscope mirror, single mirror gonioscope, mirror gonioscope, mirror gonioscope and aspheric lenses. Widely used in eye surgeries, these are precisely engineered as per specifications of the clients and provide accurate results. Features:  $\frac{1}{2}$  Affordable price  $\frac{1}{2}$  Delivery immediate  $\frac{1}{2}$  Sample Availability  $\frac{1}{2}$  High quality Specifications: Diopter  $\frac{1}{2}$ , 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50 Total 23 pieces including one red filter lens Size 37 $\frac{1}{2}$  37mm Package Aluminum Case Material Acrylic Prism Bar Set, ENT Equipment exporters, ENT Equipment manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: SITE,