



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
LENSOMETER LM- 510

Product Code :
NMCBS139030



Description :

LENSOMETER LM- 510 Crossii nes with corona target External reading with tilling facility +/- to 25 Diopters power reading 0 to 5 prism reading 0 to 1/2 to 180" cylindrical axis Batteryandmainsoperated LENSOMETER LM- 510, ENT Equipment exporters, ENT Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,