

**Product Name :**  
COMPUTERIZED LENSOMETER**Product Code :**  
NMCBS139031**Description :**

COMPUTERIZED LENSOMETER Technologically based on and supported by Zhejiang University, the company has developed and had independent intellectual property rights over full automatic Computerized Lensometer in the same step with the international mainstream. This lensometer's accuracy degree is 1 degree, so it is usually called 1 degree meter. Its technical performance has completely reached or surpasses the level of imported like products: Important Information: The only lensometer equipment at present that can measure accurately the low-degree lens below  $\pm 0.12$  in the country The lensometer equipment that can measure accurately the colored lens with a degree of transparency of 10% The lensometer equipment that can accurately and quickly measure the progressive lens. The only lensometer equipment that can display the pupil diameter Besides, the equipment can measure the pupil distance, and is equipped with UV meter and built-in thermosensitive printer It is characterized with fast measuring speed, clear interface, convenient operation, accuracy and good repeatability etc. COMPUTERIZED LENSOMETER, ENT Equipment exporters, ENT Equipment manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)