

**Product Name :**  
VIEW TESTER**Product Code :**  
NMCBS139038**Description :**

VIEW TESTER We are engaged in offering View Tester. The View Tester is distinguished by its high accuracy, wide measuring range, simplicity of operation & above all it is a really portable instrument operating on three pencil cells of size. we provide high quality tester to our clients. Features: Contour design elegant. With Sphere fast rotating disk, enabling the rapid adjustment of 3.00m-1. The unique Cross Cylinder design, convenient and flexible. Can be easily used to measure near visual acuity. VIEW TESTER, ENT Equipment exporters, ENT Equipment manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)