

Product Name :
Hospital Trial Lens Set

Product Code :
OPHTH2010001



Description :

Hospital Trial Lens Set Tender Suppliers Africa and Hospital Trial Lens Set Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Hospital Trial Lens Set Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Hospital Trial Lens Set China, Hospital Trial Lens Set Kenya.

Hospital Trial Lens Set Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Hospital Trial Lens Set Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),

Suriname (Paramaribo)

Technical Specification :

Trial frames are an ophthalmic computing instrument used by hospital ophthalmologists and opticians to check people's refractive conditions: hyperopia, nearsightedness, presbyopia, astigmatism, or color blindness.

Metal Rim: This optical lens set features metal rim, which is suitable for MRI or environments that require non-magnetic materials/equipment.

Large Aperture: With a 36mm transparent aperture, the axis marks of the cylinder lens and prism are engraved directly on the lens surface, and the axial position does not change.

Widen Slot Design: The trial lens storage slot is widened for easy removal or insertion of lenses. The standard optical lens slot effectively protects the lens.

Aluminum Case: The trial lens case is manufactured by high-quality aluminum, hard texture, combination lock, easy to carry and worry-free.

Hospital Trial Lens Set - Trial Lens Set 232 piece reduced aperture Trial Lens Set comes in a Leather case with a removable Lens tray, Lenses have Plastic Rims with Powers marked on one side. Hospital Trial Lens Set, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,