



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
RETINOSCOPE

Product Code :
NMCBS139045



Description :

RETINOSCOPE Separate control for Beam adjustment & Streak rotation 360° Beam adjustments to Divergent, Parallel & Convergent. Large Semi-Reflector interference Mirror. With one spare Bulb. Available in C & D Battery Handle, Optional Electric Handle. Streak Retinoscope in Case Streak Retinoscope Bulb Streak Retinoscope Head Streak Retinoscope Handle Streak Retinoscope with Helogen Bulb Spare Halogen Bulb. RETINOSCOPE, ENT Equipment exporters, ENT Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: MedicalHospitalEquipmentsManufacturers.com