

**Product Name :**  
Electric Obstetric Bed - Multi Function**Product Code :**  
NMCBS120013**Description :**

Electric Obstetric Bed - Multi Function This table is designed for women's obstetric and gynaecological operation, Easy operation, safe and reliable, economical and quick. Operating positions - up / down, adjustable back rest & reverse trendelenburg position. It adopts linear bearing to control the telescopic working table. Size: 2000/1450 ÷ 1/2 900/660 ÷ 1/2 620/920 mm. All motion controlled with T-Motion Power Device (Actuator) & Battery backup. Adjustable leg support and sliding grip holder. Main framework made of Epoxy Coated steel. High quality beautiful ABS engineering plastic head panel. Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward. 125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. Equipped with 8cm thick high quality foam mattress