



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Examination Couch 3 Section

Product Code :
NMCBS119007



Description :

Examination Couch 3 Section Height adjustment by electromotive drive $\pm 1/2$ The electric motor is controlled by 1 foot switch (possible option: Hand switch) $\pm 1/2$ Head part step less up- (+50 $\pm 1/2$) and downwards (-45 $\pm 1/2$) adjustable by gas spring. $\pm 1/2$ Length head part: 450 mm $\pm 1/2$ Middle part fixed. Length: 1.290 mm $\pm 1/2$ Foot part step less adjustable from downwards vertical position (-90 $\pm 1/2$) to upwards position (+28 $\pm 1/2$) by gas spring. $\pm 1/2$ Length foot part: 450 mm $\pm 1/2$ Upholstery base: Laminated wood. $\pm 1/2$ Length min.: 1.790 mm, Length max. 2.190 mm Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Piece Port of Dispatch: New Delhi Production Capacity: 100 units per month Delivery Time: As per buyer Packaging Details: Standard export packing



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,