



Email : sales@naugramedical.com

Phone: 9896600003

**Product Name :**  
Examination Table Electric

**Product Code :**  
NMCBS119008



**Description :**

Examination Table Electric Height adjustable from 560-1000mm  $\frac{1}{2}$  Scissor action provides vertical lift  $\frac{1}{2}$  Mains powered electric actuator (240 volt electrical supply)  $\frac{1}{2}$  Hand switch operation  $\frac{1}{2}$  Head section offers angle adjustment by a pneumatic pump  $\frac{1}{2}$  Paper roll holder fitted to foot end  $\frac{1}{2}$  125mm swivel braking castors  $\frac{1}{2}$  Maximum load 190kg Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Piece Port of Dispatch: New Delhi Production Capacity: 500 units per month Delivery Time: As per buyer Packaging Details: Standard export packing



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)