

Product Name :
Blood Donor Couch Electric**Product Code :**
NMCBS122004**Description :**

Blood Donor Couch Electric Motorized Donor Chairs suitable to use for blood banks, infusion treatment centers, and hospital, clinics, and aphaeresis, plasma centers and other similar discipline. Remote switch for easy repositioning All hardwood frame construction Push pillow back for donor comfort Lockable castors for easy mobility Fire Resistant Upholstery & Punishment proof laminates It has solid surface seat to ensure easy & hygienic premium rated foam, cushion and spring component for donor comfort. With the help of electric power, we can fully adjust the donor position in between head low-feet high to head high-feet lower position and arm rest (left and right) swing out for easy seating as well as up & down of chair such position minimization of donor. High quality contract grade vinyl upholstery with antimicrobial/antibacterial coating to prevent from stain, abrasion and bacterial growth. All exposed metal parts are chrome plate for easy cleaning & durability. Technical Specification: Dimension: minimum 60 inch (L) x36 inch (W) Head height reachable up to approx. 45 inch Foot height reachable up to approx. 35 inch Seat Height approx. 25 inch Capacity of patient weight. 150kg Power: 220VAC, 50Hz, Single Phase. Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Pieces Port of Dispatch: New Delhi Production Capacity: 100 units per month Delivery Time: As per buyer Packaging Details: Standard export packing

