

**Product Name :**  
Blood Donor Couch**Product Code :**  
NMCBS122005**Description :**

Blood Donor Couch Blood sampling chair with fixed height, Height= 670 mm Comfortable upholstery with good density, width 550 Stainless Steel quality, epoxy paint Adjusting back rest via self locking mechanism system Adjusting leg rest via a self locking mechanism system Equipped with 2 blood test splints, every position adjustable by serrated blocks system Load capacity: Dynamic 150kg Other Information Pay Mode Terms: T/T (Bank Transfer) Minimum Order Quantity: 1 Piece Port of Dispatch: New Delhi Production Capacity: 500 units per month Delivery Time: As per buyer Packaging Details: Standard packing

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)