

**Product Name :**  
ICU Bed 2 function Electric**Product Code :**  
NMHOPB00017**Technical Specification :**

ICU Bed 2 function Electric with electrically operated remote controlled back rest 0-70°½, knee rest 0-35°½ position. All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Size: L 2140mm x ½ W 940mm x ½ H 500mm. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors (all with individual brakes.) Provision for I.V. Rod on both sides of the bed