

Product Name :
Naugramedical ICU Bed Manual Three Function

Product Code :
NMHOPB00018



Description :

Naugramedical ICU Bed Manual Three Function India China

Technical Specification :

Naugramedical ICU Bed Manual Three Function with manually operated back rest 0-80

knee rest 0-35

and height position. Size: L 2160mm

W 940mm

H 530-810mm (Adjustable Height). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]