

**Product Name :**  
ICU Bed Manual Back Rest 3 Function**Product Code :**  
NMHOPB00021**Description :****Naugramedical ICU Bed Manual Back Rest 3 Function India China****Technical Specification :**

Naugramedical ICU Bed Manual Three Function with manually operated back rest 0-80, knee rest 0-35 and height position.

Size: L 2140mm

W 940mm

H 530-720mm (Adjustable Height). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. (B type design) Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. Bumper protection at 4 corners. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed.

