

**Product Name :**  
ICU Bed Electric IV Rod Five Function

**Product Code :**  
NMHOPB00022



**Description :**

**ICU Bed Electric IV Rod Five Function India China**

**Technical Specification :**

I.C.U. Bed, Electric, Five Function ICU Bed (5 function) with electrically operated remote controlled back rest tilting 0-80

knee rest tilting 0-35

trendelenburg tilting 0-20

Reverse trendelenburg tilting 0-20

Size: L 2160mm

W 940mm

H 480-730mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with CPR device which allows nursing staff & doctors to simultaneously lower the back rest to a full flat position from any degree of elevation in order to save vital time in critical situation. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



## Medical Hospital Equipments Manufacturers

**Website:**Medical Hospital Equipments Manufacturers,