

Product Name :
Medical ICU Bed Electric Five Function

Product Code :
NMHOPB00020



Description :

Medical ICU Bed Electric Five Function India China

Technical Specification :

I.C.U. Bed, Electric, Five Function ICU Bed (5 function) with electrically operated remote controlled back rest tilting 0-80

knee rest tilting 0-35

trendelenburg tilting 0-20

Reverse trendelenburg tilting 0-2

Size: L 2160mm

W 940mm

H 480-730mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,