

Product Name :
ICU Bed Electric Four Function With Sitting Position

Product Code :
NMHOPB00025



Description :

ICU Bed Electric Four Function With Sitting Position India China

Technical Specification :

I.C.U. Bed, Electric, Four Function With Sitting Position ICU Bed (4 function) with electrically operated remote controlled sitting position 0-70

back rest tilting 0-70

knee rest tilting 0-35

Size: L 2140mm

W 940mm

H 570-780mm (Adjustable Height). All functions controlled with electric Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. With Phenolic Laminating Sheet for X-ray. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section mesh mattress platform top. Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,