

Product Name :
Hospital ICU Bed Electric Six Function With Sitting Position

Product Code :
NMHOPB00026



Description :

Hospital ICU Bed Electric Six Function With Sitting Position

Technical Specification :

I.C.U. Bed, Electric, Six Function With Sitting Position ICU Bed (6 function) with electrically operated sitting position 0-70

back rest tilting 0 - 75%, knee rest tilting 0-35

trendelenburg tilting 0 - 15

Reverse trendelenburg tilting 0 - 15

Size: 2155 x 950 x 570 / 780 mm (Adjustable Height). All functions controlled with Electric Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]