

**Product Name :**  
Naugramedical ICU Bed Electric Three Function

**Product Code :**  
NMHOPB00027



**Description :**

**Naugramedical ICU Bed Electric Three Function India China**

**Technical Specification :**

Naugramedical ICU Bed Electric Three Function with electrically operated remote controlled back rest 0-70, knee rest 0-35 and height position. Size: L 2160mm

W 940mm

H 520-770mm (Height Adjustable). All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Easy to operate control panels on both sides of guard rails. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed

---