



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Delivery Bed

Product Code :
NMIVP20010



Description :

Delivery Bed Size : 1800L x 550W x 800H mm. Adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Supplied with 2 knee crutches, adjustable in height and width. Framework made of mild steel epoxy coated finish with upholstered top. Supplied with S.S. bowl. Knockdown construction



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,