

Product Name :
Dressing Trolley, S.S.**Product Code :**
NMTOY10027**Description :**

Dressing Trolley, S.S. Trolley made of stainless steel. Size: 660 ĩ½ 440 ĩ½ 925 mm. 2 S.S. shelves with guard rails on 3 sides. With 1 waste bin & 1 S.S. Bowl. S.S. Dressing tray attached at the top. 1 S.S. drawer below the first shelf. 4 swivel castors (2 with brakes).

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,