



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Mayo Type Trolley Height Adjustable

**Product Code :**  
NMTOY10037



**Description :**

Mayo Type Trolley Height Adjustable Trolley with Stainless Steel framework. Used in the O.T. rooms for keeping instruments during surgery. Size: 600  $\times$  400  $\times$  850 / 1000 mm. Adjustable height from 850 - 1000 mm. 2 Swivel castors.



**Medical Hospital Equipments Manufacturers**

Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)