

**Product Name :**  
Dressing Trolley**Product Code :**  
NMTOY10043**Description :**

Dressing Trolley Size : 76L x 51W x 82H cms. Frame work made of ERW tubular pipe of 31mm dia. mounted on four 10 cms. dia castors. Two fixed type S.S. Heavy gauge shelves. Supplied complete with bucket and basin. Pre-treated and epoxy powder coated. Knockdown construction.

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,