



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Dressing Trolley - Stainless steel

Product Code :
NMTOY10044



Description :

Dressing Trolley - Stainless steel Size : 76L x 51W x 82H cms. Frame work made of Stainless steel mounted on four 10 cms. dia castors. Two fixed type S.S. heavy gauge shelves. Completely supplied with bucket and basin. Knockdown construction



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,