



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Bedside Locker std. Hospital Model

Product Code :
NMLOK10009



Description :

Bedside Locker std. Hospital Model Overall size : 40L x 40W x 82H cms. Machine pressed sheet box with door of heavy gauge. Assembled to CRC tubular legs of 2.5 cms. dia. fitted with PVC stumps. Stainless steel polished top of 22G. Pre-treated and epoxy powder coated.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,