



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Bedside Locker - Deluxe

**Product Code :**  
NMLOK10010



**Description :**

Bedside Locker - Deluxe Overall size : 40L x 40W x 82H cms. Complete locker box made of 22G CRC machine pressed sheet closed from three sides and a suitable drawer for extra storage. Stainless steel top of 22G in matt finish. Complete locker mounted on two PVC Stumps in front and two 5 cms. dia castors at the back for easy moving. Pre-treated and epoxy powder coated finish



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,