



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Stretcher Army 3 Fold

Product Code :
NMSTR10015



Description :

Stretcher Army 3 Fold Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 188x45x4cm, Folded Size (LxWxH): 46x45x15cm. Self-weight: 6 kg, Load Bearing: 159 kg. Supplied with carry bag



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,