



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Stretcher Army Single Fold, with Telescopic Lifting Handles

Product Code :
NMSTR10022



Description :

Stretcher Army Single Fold, with Telescopic Lifting Handles Light Weight, easy clean, aluminium alloy framework with telescopic lifting handles. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 220x54x15cm. Folded Size (LxWxH): 180x12x15cm. Self-weight: 8 kg, Load Bearing: 159 kg.



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,