

Product Name :
Stretcher Army Double Fold With Telescopic Lifting
Handles

Product Code :
NMSTR10023



Description :

Stretcher Army Double Fold With Telescopic Lifting Handles Light Weight, easy clean, aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 220x54x15cm. Folded Size (LxWxH): 95x20x10cm. Self-weight: 8 kg, Load Bearing: 159 kg. Supplied with carry bag

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)