



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Salvage Stretcher Yellow

**Product Code :**  
NMSTR10026



**Description :**

Salvage Stretcher Yellow Designed for fractured and immobile serious patients. Light Weight, easy clean, anodized aluminium alloy framework with safety belts. Size (LxWxH): Max. unfolded Size (LxWxH): 214x42x7cm, Min. unfolded Size (LxWxH): 167x42x7cm. Detachable and folding. This Salvage stretcher can be used with head immobilizer. Length of stretcher adjustable according to patient's body. The foot-end is designed as a narrow-frame structure. Self-weight: 8 kg, Load Bearing: 159 kg



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,