

Product Name :
Head Immobilizer**Product Code :**
NMSTR10029**Description :**

Head Immobilizer Made of high-density plastic material, large ear holes for monitoring the patient's ear canal. It is used with Salvage stretcher & Spine Stretchers. Minimal interference with X-ray, MRI or CT scanning procedures. Waterproof plastic coating makes it easy to clean & prevents bacterial growth on it

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: MedicalHospitalEquipmentsManufacturers.com