



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Head Immobilizer

**Product Code :**  
NMSTR10029



**Description :**

Head Immobilizer Made of high-density plastic material, large ear holes for monitoring the patient's ear canal. It is used with Salvage stretcher & Spine Stretchers. Minimal interference with X-ray, MRI or CT scanning procedures. Waterproof plastic coating makes it easy to clean & prevents bacterial growth on it



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)