

Product Name :
Baby Bassinet with Moulded Crib**Product Code :**
NMINF10003**Description :**

Baby Bassinet with Moulded Crib.

Removable moulded Plastic Crib (without joints).

Framework made of stainless steel.

1 S.S. shelf at the bottom or Plastic Tray

4 swivel castor (all with brakes).

Size: 810 X 460 X 810 mm. Baby Bassinet with Moulded Crib , Infant Bed & Cradle Exporters, Infant Bed & Cradle Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,