

**Product Name :**  
Pediatric Bed S.Steel**Product Code :**  
NMINF10005**Description :**

Semi-Fowler Bed with mechanically operated Back rest by crank mechanism.

Steel frame work with 2 section top.

S.S. side railings with Up & Down movement.

With S.S. Shelf at the bottom for keeping blankets, basin etc.

Tubular Head & Foot Bows also made of stainless steel.

Mounted on four silent swivel castor (all with brakes).

Provision for I.V. Rod.

Freight saving knock-down construction.

Size: 1840 × 800 × 580 mm. Pediatric Bed S.Steel, Infant Bed & Cradle Exporters, Infant Bed & Cradle Manufacturers