



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Instrument Trolley

Product Code :
NMTOY10003



Description :

Instrument Trolley Frame work made of ERW tubular pipe mounted on four 7.5 cms. dia. castors.

Two fixed stainless steel shelves of heavy gauge.

Pre-treated and epoxy powder coated.

Also available :

Instrument Trolley - All stainless steel (NMTOY10002A-0086)

Instrument Trolley Without railing. (NMTOY10002B-0087)

Instrument Trolley Size: 68L x 45W x 82H cms.

Instrument Trolley Size: 76L x 45W x 82H cms.

Instrument Trolley Size: 91L x 45W x 82H cms.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,