

**Product Name :**  
Artificial Resuscitator Silicon**Product Code :**  
NMANI10001**Description :**

Artificial Resuscitator Silicon, with 2 Face Mask

Size:Adult, Child, Infant Artificial Resuscitator Silicon, Anaesthesia Equipments Exporters, Anaesthesia Equipments Manufacturers

**Technical Specification :**

Hand operated, compressible self-refilling ventilation bag with a non-rebreathing patient valve with pressure limitation valve to minimize the risk of over-inflating, intake valve with a nipple for oxygen tubing.

Supplied with oxygen reservoir bag (suitable capacity) and 2-meter oxygen tubing.

Resuscitator can be disassembled for cleaning, disinfecting or autoclaving process.