



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size)

**Product Code :**  
NMCBS35089



**Description :**

ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR is a simple, non-specific screening test utilized for monitoring the presence of inflammation in the body. ESR tubes contain tri-sodium citrate solution and are used for blood sedimentation analysis. Mixing recommendation: ESR tube should be gently inverted 180° back and forth 8-10 times



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)