



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Multi - Sample Blood Collection Needles With Hub (MOQ:
100,000 Pcs Assorted Size)

Product Code :
NMCBS35091



Description :

Multi - Sample Blood Collection Needles With Hub (MOQ: 100,000 Pcs Assorted Size) The NET Multi -Sample Blood Collection Needles With Hub with optical venipuncture control can make an important contribution to simplifying blood collection routine for phlebotomists. The needle has a transparent plastic hub, which functions as a view window. If venipuncture is successful, the view window is coloured red due to the blood flow



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,