

Product Name :
Hospital Bedside Locker DX**Product Code :**
NMLOK10002**Description :**

Complete locker box made of 22G CRC machine pressed sheet closed from three sides and a suitable drawer for extra storage.

Stainless steel top of 22G in matt finish.

Complete locker mounted on two PVC Stumps in front and two 5 cms. dia castors at the back for easy moving.

Pre-treated and epoxy powder coated finish.

Size(in) : 16 L x 16 W x 32 H Hospital Bedside Locker DX, Bedside Lockers Exporters, Bedside Lockers Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers****Website:**Medical Hospital Equipments Manufacturers,